## **General protocol: Chromogenic Immunohistochemistry (IHC)** DAY TWO

1. Prepare three 250mL wash buckets by adding **25mL of 10X PBS** to each bucket, then

topping each bucket up to the **250mL line with distilled water**.

- Add 2.5mL Tween (1%) as well to the middle bucket in the series of three washes.
- Take slides out of the humidity chamber and check to make sure none of the slides became dry overnight.
  - If this occurred, make a note of which samples this happened to and consider removing from analysis.
- 3. Move slides to a slide rack while holding the rack above the first wash bucket so that the liquid on the slides is caught by the bucket. Place in this bucket for 5 minutes.
- 4. Wash in second (Tween-containing) and then third bucket for 5 minutes each.
- 5. Set a timer for 30 minutes. Blot the first slide and apply minimum necessary amount of drops of **species-appropriate secondary antibody** onto it. Start the timer and apply secondary antibody to each subsequent slide in **regular intervals** (30 seconds, 1 minute, whatever works) keeping the slides in a line in the order they received secondary antibody.
  - Incubate at room temperature.
  - Disperse antibodies with pipette tip if necessary.
  - While secondary antibodies incubate, prepare all materials needed to finish slide prep:

- Make DAB solution by calculating the volume needed to apply 200µL
  per slide, then adding 30µL DAB enzyme for every 1000µL of DAB
  diluent needed. Invert tube repeatedly to mix
- Fill a 100mL bucket with tap water (NOT distilled this time) and place slide rack into tap water so it is ready.
- Place kimwipes or white piece of paper near this bucket for checking DAB development.
- Empty a wash bucket, rinse, and fill it to the 250mL line with tap water.
  Place beside hematoxylin in fume hood.
- Fill 2 buckets each with 95% ethanol, 100% ethanol, and xylenes.
- 6. Transfer slides into rack and into the first wash bucket at the same time interval as you applied secondary antibody.
  - This ensures that each slide incubates with secondary antibody for a full 30 minutes.
- When each slide has been placed in the rack in the first wash bucket, perform all 3 washes for 5 minutes each.
- Reset timer. Choose which slide should receive DAB first; if optimizing, choose the most concentrated antibody prep. Do not use control as first slide. Add 200µL DAB solution to this slide and start the timer.
  - As with secondary antibody application, add DAB solution to each subsequent slide at a **regular time interval**

- Check DAB development as you go! Amount of time needed with DAB is variable but usually within **2-7 minutes**.
- 9. When DAB development is sufficient, place slides into rack in the tap water bucket. Tap water stops the DAB reaction.
  - Stop DAB reaction after each slide has incubated for the same amount of time by adding them to the bucket at the same time interval DAB was applied.
  - Consider working with a partner if processing large numbers of slides.
- 10. When all slides have been placed into the tap water bucket, transfer loaded rack into the hematoxylin bucket in the fume hood and incubate for two minutes.
- 11. Move slides into bucket of water. Bring this bucket to the sink and replace water + wash slides until water is clear.
  - Do not let slides come into contact with the running tap.
- 12. Dip rack of slides into bucket of bluing buffer 7 times.
- 13. Clear slides by dipping them 15 times each in 2 buckets of 95% ethanol, 2 buckets of 100% ethanol, and 2 buckets of xylenes in that order. Leave slide rack in final xylene bucket while mounting slides
- 14. With small pestle put **Permount** on the slide encircling the sample. Squish slide cover on and press to disperse Permount on top of sample and to **push any bubbles away.** 
  - Keep in the fume hood
- 15. Wait **at least 2 hours**, then gently clean slides with kimwipe dipped in xylene.

Cox Systems Biology – coxsystemsbiology.medicine.utoronto.ca

16. Wait a few hours more until all xylenes have evaporated and mounting media is

completely hardened before viewing/imaging