

General protocol: Fluorescent Immunohistochemistry (IHC):

DAY ONE:

1. Label all slides w/ **pencil**.
 - For organoids: do not use slide warmer as it disrupts organoid integrity.
 - For tissues: slide warmer can be used but may cause sections to fall from slides more easily.
2. Dewax slides by incubate slides in series of 2 buckets of **xlenes** for 10 minutes each, 2 buckets of **100% ethanol** for 5 minutes each, and 2 buckets of **95% ethanol** for 5 minutes each.
3. Incubate slides in bucket of **distilled water** for 10 minutes
4. Move slides to a bucket of appropriate **antigen retrieval solution***, place bucket in steamer, and steam for 20 minutes to 1hr.
5. After antigen retrieval steam, allow slides to cool in freezer, then transfer to a bucket of room temperature distilled water.
 - Monitor slides/bucket during cooldown for freezing/ice formation.
 - While slides cool, dilute antibodies in diluent of **casein buffer**.
 - Casein blocks nonspecific antibody binding.
 - Aim for **200µL of diluted antibody per slide/sample**.
 - Place moistened paper towels inside of humidity chamber and flatten them within grooves to accommodate slides.
6. Fill 3 buckets w/ **250mL 1x PBS** (dilute from 10x stock). Add 1% Tween to the first bucket.
Transfer slides in rack & wash for 5 minutes in each bucket (PBS-Tween → PBS → PBS)

- *If using TLAQ:* detergents will wash away TLAQ and therefore are incompatible. All steps involving detergents should be completed prior to TLAQ application.
- *If using TLAQ:* during washes, prepare **200µL** of TLAQ working solution per slide/sample by **diluting stock solution 1:20 in 70% EtOH**.

7. Working one slide at a time:

1. Blot water off of slides at a site far from the sample.
2. Draw a wide circle (at least 1cm gap between circle + organoids) around sample with hydrophobic pen.
 - a. Do not press too hard with pen.
 - b. Ensure each circle is complete.
3. Allow pen to dry very briefly then transfer slide into humidity chamber.
4. *If using TLAQ:*
 - a. Starting with 1 slide: set timer to keep track of intervals between each slide; prepare **250mL of PBS** for wash.
 - b. Add **200µL of TLAQ** to slide.
 - i. If necessary, gently disperse antibody solution with pipette tip just contacting the solution itself, not touching the sample/slide!
 - ii. TLAQ is light sensitive. Slides should be kept from light sources as much as possible once TLAQ has been applied.

- c. Repeat for all slides while being mindful of the interval between each slide.
 - i. Specific interval does not matter as long as they are consistent between slides.
 - ii. Keep the slides in a line in the order they TLAQ for ease of tracking.
 - d. Following incubation of 3-5 minutes w/ TLAQ, transfer each slide into PBS per your time intervals and order.
 - i. This ensures that each slide incubates with secondary antibody for a full 1 hour.
 - ii. Prepare 2 additional buckets of **250mL PBS**.
 - e. Wash for 5 minutes in each of the 3 PBS buckets.
 - f. Repeat step 7.1 for each slide prior to step 7.5
5. Apply primary antibodies:
- a. Add **200µL of diluted primary antibodies** per slide/sample.
8. Place lid on humidity chamber and carefully move the chamber to the +4°C fridge for overnight incubation of primary antibodies.