

Flow Cytometric Analysis of Murine Decidual Immune Cells (E6.5-12.5)

Staining

Required Materials:

- ☐ 1.5mL Eppendorf tubes
- ☐ Zombie-series amine-reactive viability stain
 - Alternatives: DAPI, hoerscht, propidium iodide, etc.
- ☐ Antibodies
- ☐ FcX (anti-CD16/32) (Biolegend)
- ☐ Compensation beads
 - For viability: ArcReactive +/- compensation beads (Invitrogen)
 - For antibodies: UltraComp eBeads (Invitrogen)
- ☐ PBS/HBSS
- ☐ 1% w/v PBS-BSA solution
- ☐ Refrigerated centrifuge
- ☐ 4% paraformaldehyde (PFA)